

Measuring rod for children and adults

seca 216

- + Flexible measuring range options
- + Easy-to-read measurement scale
- + Adjustable head positioner
- + Choice of two units of measurement
- + Simple assembly



seca 216

Measuring rod for children and adults



Easy to read

The result can be read easily during the measurement, because the scale runs along the side of the measuring rod. A large, easily visible display window helps with the accuracy of the reading. Even after taking the measurement, the result remains constant: The head positioner can be locked with one handle and remains exactly in its position when the patient leaves his or her place on the measuring rod.

Flexible measuring range options

Whether big or small – the seca 216 is equally suitable for measuring children and adults. Because: The measuring rod has a removable scale from 1.4 to 90 inches. The scale can be flexibly inserted into a 55.3 inches long plastic profile. The user decides which measuring range appears in the profile.



Easy installation

The seca 216 can be mounted on the wall at any height. The scale itself serves as an assembly aid: If the zero point of the scale is placed on the floor and the scale is tightly fixed to the measuring rod, the desired mounting height can be easily determined. And the measured values are guaranteed to be accurate to the millimeter.

Measuring range	max. 55.3 in / 138 cm selectable from 1.4 – 90 in / 3.5 – 230 cm
Graduation	1/16 in / 1 mm
Dimensions (WxHxD)	4.7 x 59.1 x 8.5 in / 120 x 1,500 x 217 mm
Net weight	3.7 lbs / 1.7 kg

seca north america | west
13601 Benson Avenue
Chino, CA 91710
USA

seca north america | east
1439 Stuart Engals Blvd, Suite 200
Mount Pleasant, SC 29464
USA

T +1 800 542 7322
F +1 888 705 7397
E info.us@seca.com

Precision for health
seca.com